

## SUV39H1 Rabbit Polyclonal Antibody

### Catalog #: EAB11607

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Polyclonal	WB, ELISA	48	Human, Mouse, Rat

### Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB(Western Blotting)	1:500-2000
ELISA(Enzyme-linked Immunosorbent Assay)	1:5000-20000

### Product Information

Conjugate	Unconjugate
Specificity	SUV39H1 Rabbit Polyclonal Antibody detects endogenous levels of SUV39H1 protein.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks
UniProt ID	<a href="#">O43463</a>
Entrez-Gene ID	<a href="#">6839</a>

### Product Description

SUV39H1, suppressor of variegation 3-9, selectively methylates Lysine 9 of the amino terminus of Histone H3 to generate a binding site for HP1 proteins. These HP1 proteins belong to a family of heterochromatic adaptor molecules that are implicated in both gene silencing and supra-nucleosomal chromatin structure. SUV39H1 contains both SET and chromo domains and is ubiquitously expressed. The enrichment of SUV39H1 at heterochromatic foci during interphase and centromere-specific localization during metaphase depends on the C-terminal SET domain. SUV39H1 is phosphorylated specifically at the G1/S cell cycle transition and, when forcibly expressed, suppresses cell growth. SUV39H1 acts as a long-range repressor that is capable of acting over several kilobases to silence basal promoters.

**For Reserch Use Only. Not For Use In Diagnostic Procedures**