

Phospho-CaSR (Thr888) Rabbit Polyclonal Antibody

Catalog #: EAB11635

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Polyclonal	WB, IF/ICC, ELISA	121	Human, Mouse, Rat

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB (Western Blotting)	1:500-2000
IF/ICC (Immunofluorescence/Immunocytochemistry)	1:50-300
ELISA (Enzyme-linked Immunosorbent Assay)	1:5000-20000

Product Information

Conjugate	Unconjugate
Specificity	Phospho-CaSR (Thr888) Rabbit Polyclonal Antibody detects endogenous levels of CaSR only when phosphorylated at Thr888.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks
UniProt ID	P41180
Entrez-Gene ID	846

Product Description

Extracellular calcium-sensing receptor (CaSR), also designated parathyroid cell calcium-sensing receptor, is an integral membrane protein that belongs to the G protein-coupled receptor 3 family. CaSR is involved in maintaining a stable calcium concentration by acting as a sensor of the extracellular calcium levels for the parathyroid and kidney. Its activity is mediated by a G protein which activates a phosphatidylinositol-calcium second messenger system. Defects that activate CaSR cause autosomal dominant hypocalcemia, whereas mutations that inactivate the protein cause familial hypocalciuric hypercalcemia. CaSR is expressed mainly in kidney, and is also expressed in intestine, placenta and brain.

For Reserch Use Only. Not For Use In Diagnostic Procedures