

CRCP Rabbit Polyclonal Antibody

Catalog #: EAB11728

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Polyclonal	WB, ELISA	17	Human, Mouse, Rat

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB (Western Blotting)	1:500-2000
ELISA (Enzyme-linked Immunosorbent Assay)	1:5000-20000

Product Information

Conjugate	Unconjugate
Specificity	CRCP Rabbit Polyclonal Antibody detects endogenous levels of CRCP protein.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks
UniProt ID	Q75575
Entrez-Gene ID	27297

Product Description

The CGRP-receptor complement protein (CGRP-RCP) also designated CRCP is an accessory protein that is expressed in CGRP responsive tissues. CRCP is associated with the acrosome, thereby playing a role in reproduction. Specifically, CRCP couples the CGRP receptor to downstream effectors. This gene encodes a membrane protein that functions as part of a receptor complex for a small neuropeptide that increases intracellular cAMP levels.

For Reserch Use Only. Not For Use In Diagnostic Procedures