

## Phospho-MKK3 (Ser189) Rabbit Polyclonal Antibody

### Catalog #: EAB11919

| Host/Isotype | Clonality  | Applications         | MW (kDa) | Reactivity        |
|--------------|------------|----------------------|----------|-------------------|
| Rabbit IgG   | Polyclonal | WB, IHC-P, IF, ELISA | 39       | Human, Mouse, Rat |

### Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

|                                                  |              |
|--------------------------------------------------|--------------|
| <b>WB</b> (Western Blotting)                     | 1:500-2000   |
| <b>IHC-P</b> (Immunohistochemistry-Paraffin)     | 1:50-300     |
| <b>IF</b> (Immunofluorescence)                   | 1:50-300     |
| <b>ELISA</b> (Enzyme-linked Immunosorbent Assay) | 1:5000-20000 |

### Product Information

|                       |                                                                                                                                        |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Conjugate</b>      | Unconjugate                                                                                                                            |
| <b>Specificity</b>    | Phospho-MKK3 (Ser189) Rabbit Polyclonal Antibody detects endogenous levels of MKK3 only when phosphorylated at Ser189.                 |
| <b>Purification</b>   | Affinity purification                                                                                                                  |
| <b>Concentration</b>  | 1mg/ml                                                                                                                                 |
| <b>Format</b>         | Liquid                                                                                                                                 |
| <b>Formulation</b>    | In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol                                                               |
| <b>Shipping</b>       | Gel Pack                                                                                                                               |
| <b>Storage</b>        | Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks |
| <b>UniProt ID</b>     | <a href="#">P46734</a>                                                                                                                 |
| <b>Entrez-Gene ID</b> | <a href="#">5606</a>                                                                                                                   |

### Product Description

The protein encoded by this gene is a dual specificity protein kinase that belongs to the MAP kinase kinase family. This kinase is activated by mitogenic and environmental stress, and participates in the MAP kinase-mediated signaling cascade. It phosphorylates and thus activates MAPK14/p38-MAPK. This kinase can be activated by insulin, and is necessary for the expression of glucose transporter. Expression of RAS oncogene is found to result in the accumulation of the active form of this kinase, which thus leads to the constitutive activation of MAPK14, and confers oncogenic transformation of primary cells. The inhibition of this kinase is involved in the pathogenesis of Yersinia pseudotuberculosis. Multiple alternatively spliced transcript variants that encode distinct isoforms have been reported for this gene.

**For Reserch Use Only. Not For Use In Diagnostic Procedures**