

SGK3 Rabbit Polyclonal Antibody

Catalog #: EAB13344

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Polyclonal	WB, ELISA	57	Human, Mouse, Rat

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB(Western Blotting)	1:500-2000
ELISA(Enzyme-linked Immunosorbent Assay)	1:5000-20000

Product Information

Conjugate	Unconjugate
Specificity	SGK3 Rabbit Polyclonal Antibody detects endogenous levels of SGK3 protein.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol.
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks.
UniProt ID	Q96BR1
Entrez-Gene ID	23678

Product Description

Serine/threonine-protein kinase Sgk3 (SGK3), also designated serum/glucocorticoid regulated kinase 3, belongs to the Ser/Thr protein kinase family of proteins. The serum- and glucocorticoid-regulated kinase proteins are closely related to the Akt protein family. SGK1, a homolog of SGK3, activates ion channels, in particular potassium (K⁺) channels. SGK2 and SGK3 have been found to also be involved in this activation process, making all three of these proteins important regulators for cell proliferation, epithelial transport and neuromuscular excitability. SGK3 acts as a mediator of IL-3 dependent survival signals in the cell. It localizes to the early endosome and in vesicle-like structures. SGK3 is a widely expressed protein, but it is primarily detected in kidney, liver, pancreas, brain and heart. Phosphorylation of SGK3 at residue Ser 486 leads to an increase in SGK3 activation.

For Reserch Use Only. Not For Use In Diagnostic Procedures