

Lamin B1 Rabbit Polyclonal Antibody

Catalog #: EAB13626

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Polyclonal	WB, IHC-P, IF/ICC	68	Human, Mouse, Rat

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB(Western Blotting)	1:500-2000
IHC-P(Immunohistochemistry-Paraffin)	1:50-300
IF/ICC(Immunofluorescence/Immunocytochemistry)	1:50-300

Product Information

Conjugate	Unconjugate
Specificity	Lamin B1 Rabbit Polyclonal Antibody detects endogenous levels of Lamin B1 protein.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks
UniProt ID	P20700
Entrez-Gene Id	4001

Product Description

Lamin B1 is a 586 amino acid protein that is encoded by a gene which, when mutated, is involved in the pathogenesis of autosomal dominant adult-onset leukodystrophy (ADLD), a disease characterized by cerebellar dysfunction and symmetric demyelination of the central nervous system. B-type lamins undergo a series of modifications, such as farnesylation and phosphorylation. Increased phosphorylation of the lamins occurs before envelope disintegration and probably plays a role in regulating lamin associations. Nuclear Lamin B is fragmented as a consequence of apoptosis by an unidentified member of the ICE family.

For Reserch Use Only. Not For Use In Diagnostic Procedures