

Clathrin heavy chain/CLTC Rabbit Polyclonal Antibody

Catalog #: EAB13761

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Polyclonal	WB, IF/ICC	192	Human, Mouse, Rat

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB(Western Blotting)	1:500-2000
IF/ICC(Immunofluorescence/Immunocytochemistry)	1:50-300

Product Information

Conjugate	Unconjugate
Specificity	Clathrin heavy chain/CLTC Rabbit Polyclonal Antibody detects endogenous levels of Clathrin heavy chain/CLTC protein.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks
UniProt ID	Q00610
Entrez-Gene ID	1213

Product Description

Clathrin is a major cytosolic coat protein in pits and vesicles originating from the plasma membrane and the trans-Golgi network. In receptor-mediated endocytosis, receptor proteins are engulfed by clathrin-coated vesicles. Clathrin is composed of three heavy chains and three light chains which associate non-covalently to form a triskelion structure. Clathrin heavy chain (HC) is composed of a terminal globular domain, a distal segment and a proximal segment containing a light chain binding site. The proximal segment of the Clathrin HC protein is essential for interactions between clathrin heavy chains and light chains which result in the formation of the triskelion structure.

For Reserch Use Only. Not For Use In Diagnostic Procedures