

DNAJC15 Rabbit Monoclonal Antibody

Catalog #: EAB21468

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Monoclonal	WB, IHC-P	17	Human, Mouse, Rat

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB(Western Blotting)	1:500-2000
IHC-P(Immunohistochemistry-Paraffin)	1:50-200

Product Information

Conjugate	Unconjugate
Specificity	DNAJC15 Rabbit Monoclonal Antibody detects endogenous levels of DNAJC15 protein.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks
UniProt ID	Q9Y5T4
Entrez-Gene Id	29103

Product Description

The DnaJ family is one of the largest of all the chaperone families and has evolved with diverse cellular localization and functions. The presence of a J domain defines a protein as a member of the DnaJ family. DnaJ heat shock induced proteins are from the bacterium *Escherichia coli* and are under the control of the htpR regulatory protein. The DnaJ proteins play a critical role in the HSP 70 chaperone machine by interacting with HSP 70 to stimulate ATP hydrolysis. DnaJ proteins are important mediators of proteolysis and are involved in the regulation of protein degradation, exocytosis and endocytosis. MCJ (methylation-controlled J protein), also known as HSD18, DNAJD1 or DNAJC15, is a 150 amino acid ubiquitously expressed single-pass membrane protein containing one J domain. Localizing to the golgi apparatus and only present in vertebrates, MCJ may be associated with increased chemotherapeutic resistance in ovarian cancer by inducing expression of the Mdr drug transporter and preventing intracellular drug accumulation.

For Reserch Use Only. Not For Use In Diagnostic Procedures